


New version of Segment CMR

Today we are releasing version 2.1 of Segment CMR with reworked, cleaner and improved user interface. The new interface is based on the result of a usability study we recently made.



List of main improvements

- LV and RV Volumes are now always displayed in the report panel and not dependent on active image stack.
- Larger and clearer icons. All icons are moved to one location.
- New row of icons that depend on active mode
- More functions that are available from icons
- Flow graph is displayed in report panel
- Import contours are automatically aligned (import cine contours to LGE for example)
- Additional hot keys.
- much more ...

Have a personal meeting with us.

We recommend all sites to have a personal meeting with us, where we can show the latest features, answer questions and help you to improve your workflow. Meetings are typically help using Webex or Skype. Please contact [Einar Heiberg](mailto:Einar.Heiberg@medviso.com) to book your meeting.

Segment CMR Newsletter

You are receiving this newsletter since you are registered as a user of the clinical image analysis software Segment CMR provided by Medviso AB. If you do not want to receive further newsletters please use the link below to unsubscribe. If you have a colleague who also would like to get this newsletter, forward it and notify me by reply to this email.

Best regards

Einar Heiberg, Founder Medviso AB
Einar@medviso.com